



ATG14 blocking peptide (CDBP5121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ATG14 autophagy related 14 [Homo sapiens (human)]
Official Symbol	ATG14
Synonyms	ATG14; autophagy related 14; ATG14L; BARKOR; KIAA0831; beclin 1-associated autophagy-related key regulator; Beclin 1-Interacting protein; ATG14 autophagy related 14 homolog; autophagy-related protein 14-like protein
Entrez Gene ID	22863
mRNA Refseq	NM_014924
Protein Refseq	NP_055739
UniProt ID	Q6ZNE5

Pathway	Regulation of autophagy; Senescence and Autophagy
Function	protein binding