



ATG2B blocking peptide (CDBP5123)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ATG2B autophagy related 2B [Homo sapiens (human)]
Official Symbol	ATG2B
Synonyms	ATG2B; autophagy related 2B; C14orf103; autophagy-related protein 2 homolog B; ATG2 autophagy related 2 homolog B
Entrez Gene ID	55102
mRNA Refseq	NM_018036
Protein Refseq	NP_060506
UniProt ID	Q96BY7