



ATOH8 blocking peptide (CDBP5128)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ATOH8 atonal homolog 8 (Drosophila) [Homo sapiens (human)]
Official Symbol	ATOH8
Synonyms	ATOH8; atonal homolog 8 (Drosophila); HATH6; bHLHa21; protein atonal homolog 8; atonal homolog 6; helix-loop-helix protein hATH-6; class A basic helix-loop-helix protein 21; basic helix loop helix transcription factor 6; basic helix-loop-helix transcription factor 6
Entrez Gene ID	84913
mRNA Refseq	NM_032827
Protein Refseq	NP_116216
UniProt ID	Q96SQ7

FunctionDNA binding; protein dimerization activity
