



BCAP29 blocking peptide (CDBP5161)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BCAP29 B-cell receptor-associated protein 29 [Homo sapiens (human)]
Official Symbol	BCAP29
Synonyms	BCAP29; B-cell receptor-associated protein 29; BAP29; BCR-associated protein 29
Entrez Gene ID	55973
mRNA Refseq	NM_001008405
Protein Refseq	NP_001008405
UniProt ID	Q9UHQ4