



CCP110 blocking peptide (CDBP5168)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BCAS3 breast carcinoma amplified sequence 3 [Homo sapiens (human)]
Official Symbol	CCP110
Synonyms	BCAS3; breast carcinoma amplified sequence 3; MAAB; GAOB1; breast carcinoma-amplified sequence 3; protein Maab1; BCAS4/BCAS3 fusion; metastasis associated antigen of breast cancer; breast carcinoma amplified sequence 4/3 fusion protein
Entrez Gene ID	54828
mRNA Refseq	NM_001099432
Protein Refseq	NP_001092902
UniProt ID	Q9H6U6

Functionprotein binding
