



CLEC2A blocking peptide (CDBP5169)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BCAS4 breast carcinoma amplified sequence 4 [Homo sapiens (human)]
Official Symbol	CLEC2A
Synonyms	BCAS4; breast carcinoma amplified sequence 4; CNOL; breast carcinoma-amplified sequence 4; BCAS4/BCAS3 fusion; breast carcinoma amplified sequence 4/3 fusion protein
Entrez Gene ID	55653
mRNA Refseq	NM_001010974
Protein Refseq	NP_001010974
UniProt ID	Q8TDM0