



BCL9L blocking peptide (CDBP5175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BCL9L B-cell CLL/lymphoma 9-like [Homo sapiens (human)]
Official Symbol	BCL9L
Synonyms	BCL9L; B-cell CLL/lymphoma 9-like; BCL9-2; DLNB11; B-cell CLL/lymphoma 9-like protein; protein BCL9-2; BCL9-like protein; B-cell lymphoma 9-like protein; nuclear co-factor of beta-catenin signalling
Entrez Gene ID	283149
mRNA Refseq	NM_182557
Protein Refseq	NP_872363
UniProt ID	Q86UU0

Pathway	Disease; RNF mutants show enhanced WNT signaling and proliferation; Signal Transduction; Signaling by WNT in cancer; Signaling by Wnt; TCF dependent signaling in response to WNT; XAV939 inhibits tankyrase; deactivation of the beta-catenin transactivating complex
Function	beta-catenin binding; transcription coactivator activity