



BCL2L13 blocking peptide (CDBP5178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BCL2L13 BCL2-like 13 (apoptosis facilitator) [Homo sapiens (human)]
Official Symbol	BCL2L13
Synonyms	BCL2L13; BCL2-like 13 (apoptosis facilitator); MIL1; BCL-RAMBO; Bcl2-L-13; bcl-2-like protein 13
Entrez Gene ID	23786
mRNA Refseq	NM_001270726
Protein Refseq	NP_001257655

UniProt ID	Q9BXK5
Pathway	Legionellosis
Function	cysteine-type endopeptidase activator activity involved in apoptotic process