



C19orf80 blocking peptide (CDBP5183)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	C19orf80 chromosome 19 open reading frame 80 [Homo sapiens (human)]
Official Symbol	C19orf80
Synonyms	C19orf80; chromosome 19 open reading frame 80; RIFL; TD26; ANGPTL8; PRO1185; PVPA599; betatrophin; lipasin; betatrophin variant 1; betatrophin variant 2; angiopoietin-like protein 8; hepatocellular carcinoma-associated gene TD26; hepatocellular carcinoma-associated protein TD26
Entrez Gene ID	55908
mRNA Refseq	NM_018687
Protein Refseq	NP_061157
UniProt ID	Q6UXH0

Functionhormone activity; protein binding
