



PTRH2 blocking peptide (CDBP5196)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PTRH2 peptidyl-tRNA hydrolase 2 [Homo sapiens (human)]
Official Symbol	PTRH2
Synonyms	PTRH2; peptidyl-tRNA hydrolase 2; BIT1; PTH2; CFAP37; CGI-147; peptidyl-tRNA hydrolase 2, mitochondrial; PTH 2; bcl-2 inhibitor of transcription 1; cilia and flagella associated protein 37
Entrez Gene ID	51651
mRNA Refseq	NM_016077
Protein Refseq	NP_057161
UniProt ID	Q9Y3E5
Pathway	Apoptosis Modulation and Signaling

Functionaminoacyl-tRNA hydrolase activity; protein binding
