



VPS13B blocking peptide (CDBP5208)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	HAPLN2 hyaluronan and proteoglycan link protein 2 [Homo sapiens (human)]
Official Symbol	VPS13B
Synonyms	HAPLN2; hyaluronan and proteoglycan link protein 2; BRAL1; brain link protein 1; brain link protein-1
Entrez Gene ID	60484
mRNA Refseq	NM_021817
Protein Refseq	NP_068589
UniProt ID	Q9GZV7
Function	hyaluronic acid binding