



BRE blocking peptide (CDBP5210)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BRE brain and reproductive organ-expressed (TNFRSF1A modulator) [Homo sapiens (human)]
Official Symbol	BRE
Synonyms	BRE; brain and reproductive organ-expressed (TNFRSF1A modulator); BRCC4; BRCC45; BRCA1-A complex subunit BRE; BRCA1/BRCA2-containing complex subunit 45; BRCA1/BRCA2-containing complex, subunit 4; brain and reproductive organ-expressed protein
Entrez Gene ID	9577
mRNA Refseq	NM_001261840
Protein Refseq	NP_001248769
UniProt ID	Q9NXR7

Function	peroxisome targeting sequence binding; polyubiquitin binding; protein binding; tumor necrosis factor receptor binding
-----------------	---
