



Bst2 blocking peptide (CDBP5214)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Bst2 bone marrow stromal cell antigen 2 [Mus musculus (house mouse)]
Official Symbol	Bst2
Synonyms	Bst2; bone marrow stromal cell antigen 2; GREG; Bst-2; CD317; C87040; DAMP-1; 2310015I10Rik; bone marrow stromal antigen 2; tetherin; HM1.24 antigen; DAMP-1 protein homolog
Entrez Gene ID	69550
mRNA Refseq	NM_198095
Protein Refseq	NP_932763
UniProt ID	Q8R2Q8