



C2orf40 blocking peptide (CDBP5222)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	C2orf40 chromosome 2 open reading frame 40 [Homo sapiens (human)]
Official Symbol	C2orf40
Synonyms	C2orf40; chromosome 2 open reading frame 40; ECRG4; augurin; esophageal cancer related gene 4 protein; esophageal cancer-related gene 4 protein
Entrez Gene ID	84417
mRNA Refseq	NM_032411
Protein Refseq	NP_115787
UniProt ID	Q9H1Z8