



C9orf116 blocking peptide (CDBP5223)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	C9orf116 chromosome 9 open reading frame 116 [Homo sapiens (human)]
Official Symbol	C9orf116
Synonyms	C9orf116; chromosome 9 open reading frame 116; PIERCE1; RbEST47; UPF0691 protein C9orf116; pierce 1; p53-induced expression in RB-null cells protein 1; p53-induced expression 1
Entrez Gene ID	138162
mRNA Refseq	NM_001048265
Protein Refseq	NP_001041730
UniProt ID	Q5BN46