



CALN1 blocking peptide (CDBP5225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein with high similarity to the calcium-binding proteins of the calmodulin family. The encoded protein contains two EF-hand domains and potential calcium-binding sites. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CALN1 calneuron 1 [Homo sapiens (human)]
Official Symbol	CALN1
Synonyms	CALN1; calneuron 1; CABP8; calcium-binding protein 8; calneuron I; calneuron-1; calcium-binding protein CABP8
Entrez Gene ID	83698
mRNA Refseq	NM_001017440

Protein Refseq	NP_001017440
UniProt ID	Q9BXU9
Function	calcium ion binding