



CASKIN1 blocking peptide (CDBP5241)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CASKIN1 CASK interacting protein 1 [Homo sapiens (human)]
Official Symbol	CASKIN1
Synonyms	CASKIN1; CASK interacting protein 1; ANKS5A; caskin-1; CASK-interacting protein 1
Entrez Gene ID	57524
mRNA Refseq	NM_020764
Protein Refseq	NP_065815
UniProt ID	Q8WXD9
Function	protein binding