



CASP4 blocking peptide (CDBP5249)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CASP4 caspase 4, apoptosis-related cysteine peptidase [Bos taurus (cattle)]
Official Symbol	CASP4
Synonyms	CASP4; CASP13; caspase-4; ERICE; CASP-4; CASP-13; caspase 13; caspase-13; caspase 4, apoptosis-related cysteine protease; evolutionary related interleukin-1-beta-converting enzyme
Entrez Gene ID	338039
mRNA Refseq	NM_176638
Protein Refseq	NP_788811
UniProt ID	O75601
Pathway	Apoptosis

Function

cysteine-type endopeptidase activity; cysteine-type endopeptidase activity
