



Casp12 blocking peptide (CDBP5258)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Casp12 caspase 12 [Mus musculus (house mouse)]
Official Symbol	Casp12
Synonyms	Casp12; caspase 12; caspase-12; CASP-12
Entrez Gene ID	12364
mRNA Refseq	NM_009808
Protein Refseq	NP_033938
UniProt ID	O08736
Pathway	Alzheimer's disease; Amyotrophic lateral sclerosis (ALS); Apoptosis; Apoptosis signaling pathway; Hepatitis B; Prion diseases; Protein processing in endoplasmic reticulum

Function

cysteine-type endopeptidase activity; cysteine-type peptidase activity; endopeptidase activity; hydrolase activity; peptidase activity; protease binding; protein binding
