



CCDC106 blocking peptide (CDBP5265)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CCDC106 coiled-coil domain containing 106 [Homo sapiens (human)]
Official Symbol	CCDC106
Synonyms	CCDC106; coiled-coil domain containing 106; ZNF581; HSU79303; coiled-coil domain-containing protein 106; protein predicted by clone 23882
Entrez Gene ID	29903
mRNA Refseq	NM_013301
Protein Refseq	NP_037433
UniProt ID	Q9BWC9
Function	protein binding