



CCDC134 blocking peptide (CDBP5266)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CCDC134 coiled-coil domain containing 134 [Homo sapiens (human)]
Official Symbol	CCDC134
Synonyms	CCDC134; coiled-coil domain containing 134; dJ821D11.3; coiled-coil domain-containing protein 134
Entrez Gene ID	79879
mRNA Refseq	NM_024821
Protein Refseq	NP_079097
UniProt ID	Q9H6E4