



CCDC141 blocking peptide (CDBP5267)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CCDC141 coiled-coil domain containing 141 [Homo sapiens (human)]
Official Symbol	CCDC141
Synonyms	CCDC141; coiled-coil domain containing 141; CAMDI; coiled-coil domain-containing protein 141; coiled-coil protein associated with myosin II and DISC1
Entrez Gene ID	285025
mRNA Refseq	NM_173648
Protein Refseq	NP_775919
UniProt ID	Q6ZP82
Function	protein binding