



CCDC47 blocking peptide (CDBP5269)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CCDC47 coiled-coil domain containing 47 [Homo sapiens (human)]
Official Symbol	CCDC47
Synonyms	CCDC47; coiled-coil domain containing 47; GK001; MSTP041; coiled-coil domain-containing protein 47
Entrez Gene ID	57003
mRNA Refseq	NM_020198
Protein Refseq	NP_064583
UniProt ID	Q96A33
Function	calcium ion binding; poly(A) RNA binding; protein binding