



CCDC69 blocking peptide (CDBP5271)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CCDC69 coiled-coil domain containing 69 [Homo sapiens (human)]
Official Symbol	CCDC69
Synonyms	CCDC69; coiled-coil domain containing 69; coiled-coil domain-containing protein 69
Entrez Gene ID	26112
mRNA Refseq	NM_015621
Protein Refseq	NP_056436
UniProt ID	A6NI79