



# FAM175A blocking peptide (CDBP5272)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

## GENE INFORMATION

Gene Name	<a href="#">FAM175A family with sequence similarity 175, member A [ Homo sapiens (human) ]</a>
Official Symbol	FAM175A
Synonyms	FAM175A; family with sequence similarity 175, member A; ABRA1; CCDC98; BRCA1-A complex subunit Abraxas; abraxas protein; coiled-coil domain containing 98; coiled-coil domain-containing protein 98
Entrez Gene ID	<a href="#">84142</a>
mRNA Refseq	<a href="#">NM_139076</a>
Protein Refseq	<a href="#">NP_620775</a>
UniProt ID	Q6UWZ7

<b>Function</b>	polyubiquitin binding; protein binding
-----------------	--

---