



CDNF blocking peptide (CDBP5288)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CDNF cerebral dopamine neurotrophic factor [Homo sapiens (human)]
Official Symbol	CDNF
Synonyms	CDNF; cerebral dopamine neurotrophic factor; ARMETL1; ARMET-like protein 1; conserved dopamine neurotrophic factor; arginine-rich, mutated in early stage tumors-like 1
Entrez Gene ID	441549
mRNA Refseq	NM_001029954
Protein Refseq	NP_001025125
UniProt ID	Q49AH0
Function	growth factor activity