



CENPW blocking peptide (CDBP5292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CENPW centromere protein W [Homo sapiens (human)]
Official Symbol	CENPW
Synonyms	CENPW; centromere protein W; CUG2; CENP-W; C6orf173; cancer-upregulated gene 2; cancer-up-regulated gene 2 protein
Entrez Gene ID	387103
mRNA Refseq	NM_001012507
Protein Refseq	NP_001012525
UniProt ID	Q5EE01
Function	DNA binding; protein binding; protein heterodimerization activity