



CHCHD7 blocking peptide (CDBP5294)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CHCHD7 coiled-coil-helix-coiled-coil-helix domain containing 7 [Homo sapiens (human)]
Official Symbol	CHCHD7
Synonyms	CHCHD7; coiled-coil-helix-coiled-coil-helix domain containing 7; COX23; coiled-coil-helix-coiled-coil-helix domain-containing protein 7; COX23 cytochrome c oxidase assembly homolog
Entrez Gene ID	79145
mRNA Refseq	NM_001011667
Protein Refseq	NP_001011667
UniProt ID	Q9BUK0