



## CIDEB blocking peptide (CDBP5300)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                         |
|-----------------|-----------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">CIDEB cell death-inducing DFFA-like effector b [ Homo sapiens (human) ]</a> |
| Official Symbol | CIDEB                                                                                   |
| Synonyms        | CIDEB; cell death-inducing DFFA-like effector b; cell death activator CIDE-B            |
| Entrez Gene ID  | <a href="#">27141</a>                                                                   |
| mRNA Refseq     | <a href="#">NM_014430</a>                                                               |
| Protein Refseq  | <a href="#">NP_055245</a>                                                               |
| UniProt ID      | Q9UHD4                                                                                  |
| Function        | identical protein binding; protein binding                                              |