



CIDEB blocking peptide (CDBP5302)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CIDEB cell death-inducing DFFA-like effector b [Homo sapiens (human)]
Official Symbol	CIDEB
Synonyms	CIDEB; cell death-inducing DFFA-like effector b; cell death activator CIDE-B
Entrez Gene ID	27141
mRNA Refseq	NM_014430
Protein Refseq	NP_055245
UniProt ID	Q9UHD4
Function	identical protein binding; protein binding