



UBQLN4 blocking peptide (CDBP5307)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	UBQLN4 ubiquilin 4 [Homo sapiens (human)]
Official Symbol	UBQLN4
Synonyms	UBQLN4; ubiquilin 4; A1U; A1Up; UBIN; CIP75; C1orf6; ubiquilin-4; connexin43-interacting protein of 75 kDa; ataxin-1 interacting ubiquitin-like protein; ataxin-1 ubiquitin-like interacting protein; ataxin-1 ubiquitin-like-interacting protein A1U
Entrez Gene ID	56893
mRNA Refseq	NM_020131
Protein Refseq	NP_064516
UniProt ID	Q9NRR5

Pathway	Protein processing in endoplasmic reticulum
Function	identical protein binding; polyubiquitin binding; protein binding