



CABS1 blocking peptide (CDBP5315)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CABS1 calcium-binding protein, spermatid-specific 1 [Homo sapiens (human)]
Official Symbol	CABS1
Synonyms	CABS1; calcium-binding protein, spermatid-specific 1; CLPH; C4orf35; NYD-SP26; calcium-binding and spermatid-specific protein 1; casein-like phosphoprotein; testis development protein NYD-SP26; calcium-binding protein, sperm-specific 1
Entrez Gene ID	85438
mRNA Refseq	NM_033122
Protein Refseq	NP_149113
UniProt ID	Q96KC9

Functioncalcium ion binding
