



CMYA5 blocking peptide (CDBP5317)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CMYA5 cardiomyopathy associated 5 [Homo sapiens (human)]
Official Symbol	CMYA5
Synonyms	CMYA5; cardiomyopathy associated 5; SPRYD2; TRIM76; C5orf10; cardiomyopathy-associated protein 5; myospryn; genethonin-3; 2310076E16Rik; dystrobrevin-binding protein 2; tripartite motif-containing 76; SPRY domain-containing protein 2; tripartite motif-containing protein 76
Entrez Gene ID	202333
mRNA Refseq	NM_153610
Protein Refseq	NP_705838
UniProt ID	Q8N3K9

Function	identical protein binding; protein binding
-----------------	--
