



COMMD7 blocking peptide (CDBP5321)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	COMMD7 COMM domain containing 7 [Homo sapiens (human)]
Official Symbol	COMMD7
Synonyms	COMMD7; COMM domain containing 7; C20orf92; dJ1085F17.3; COMM domain-containing protein 7
Entrez Gene ID	149951
mRNA Refseq	NM_001099339
Protein Refseq	NP_001092809
UniProt ID	Q86VX2
Function	NF-kappaB binding; protein binding