



CORO7 blocking peptide (CDBP5322)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Coronins, such as CORO7, constitute an evolutionarily conserved family of WD-repeat actin-binding proteins. CORO7 plays a role in Golgi complex morphology and function (Rybakin et al., 2004, 2006 [PubMed 15327992] [PubMed 16905771]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CORO7 coronin 7 [Homo sapiens (human)]
Official Symbol	CORO7
Synonyms	CORO7; coronin 7; CRN7; POD1; 0610011B16Rik; coronin-7; 70 kDa WD repeat tumor rejection antigen homolog
Entrez Gene ID	79585
mRNA Refseq	NM_001201472
Protein Refseq	NP_001188401

UniProt ID	P57737
Pathway	Metabolism of proteins; Mitochondrial protein import
Function	actin binding; NOT actin filament binding; protein binding