



KCP blocking peptide (CDBP5327)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KCP kielin/chordin-like protein [Homo sapiens (human)]
Official Symbol	KCP
Synonyms	KCP; kielin/chordin-like protein; KCP1; CRIM2; NET67; KCP-1; CRIM-2; kielin; cysteine-rich BMP regulator 2; kielin/chordin-like protein 1; cysteine-rich motor neuron 2 protein; cysteine rich BMP regulator 2 (chordin-like)
Entrez Gene ID	375616
mRNA Refseq	NM_001135914
Protein Refseq	NP_001129386
UniProt ID	Q6ZWJ8