



CRISP2 blocking peptide (CDBP5328)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CRISP2 cysteine-rich secretory protein 2 [Homo sapiens (human)]
Official Symbol	CRISP2
Synonyms	CRISP2; cysteine-rich secretory protein 2; CT36; TPX1; TSP1; GAPDL5; CRISP-2; cancer/testis antigen 36; testis-specific protein TPX-1; testis specific protein 1 (probe H4-1 p3-1); glyceraldehyde-3-phosphate dehydrogenase-like 5
Entrez Gene ID	7180
mRNA Refseq	NM_001142407
Protein Refseq	NP_001135879
UniProt ID	P16562