



# C1QTNF1 blocking peptide (CDBP5333)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                                                                                                                                                                                      |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">C1QTNF1 C1q and tumor necrosis factor related protein 1 [ Homo sapiens (human) ]</a>                                                                                                                                     |
| Official Symbol | C1QTNF1                                                                                                                                                                                                                              |
| Synonyms        | C1QTNF1; C1q and tumor necrosis factor related protein 1; GIP; CTRP1; ZSIG37; complement C1q tumor necrosis factor-related protein 1; G protein coupled receptor interacting protein; g protein-coupled receptor-interacting protein |
| Entrez Gene ID  | <a href="#">114897</a>                                                                                                                                                                                                               |
| mRNA Refseq     | <a href="#">NM_030968</a>                                                                                                                                                                                                            |
| Protein Refseq  | <a href="#">NP_112230</a>                                                                                                                                                                                                            |
| UniProt ID      | Q9BXJ1                                                                                                                                                                                                                               |

**Function**collagen binding; protein binding

---