



C1QTNF4 blocking peptide (CDBP5337)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	C1QTNF4 C1q and tumor necrosis factor related protein 4 [Homo sapiens (human)]
Official Symbol	C1QTNF4
Synonyms	C1QTNF4; C1q and tumor necrosis factor related protein 4; CTRP4; ZACRP4; complement C1q tumor necrosis factor-related protein 4; complement-c1q tumor necrosis factor-related protein 4
Entrez Gene ID	114900
mRNA Refseq	NM_031909
Protein Refseq	NP_114115
UniProt ID	Q9BXJ3

Function	cytokine activity
----------	-------------------
