



C1QTNF5 blocking peptide (CDBP5340)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of a family of proteins that function as components of basement membranes and may play a role in cell adhesion. Mutations in this gene have been associated with late-onset retinal degeneration. The protein may be encoded by either a bicistronic transcript including sequence from the upstream membrane frizzled-related protein gene (MFRP), or by a monocistronic transcript expressed from an internal promoter. [provided by RefSeq, Jun 2013]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	C1QTNF5 C1q and tumor necrosis factor related protein 5 [Homo sapiens (human)]
Official Symbol	C1QTNF5
Synonyms	C1QTNF5; C1q and tumor necrosis factor related protein 5; CTRP5; complement C1q tumor necrosis factor-related protein 5; myonectin; C1q TNF-alpha-related protein 5
Entrez Gene ID	114902

mRNA Refseq [NM_001278431](#)

Protein Refseq [NP_001265360](#)

UniProt ID Q9BXJ0
