



C1QTNF7 blocking peptide (CDBP5344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	C1QTNF7 C1q and tumor necrosis factor related protein 7 [Homo sapiens (human)]
Official Symbol	C1QTNF7
Synonyms	C1QTNF7; C1q and tumor necrosis factor related protein 7; CTRP7; ZACRP7; complement C1q tumor necrosis factor-related protein 7; complement-c1q tumor necrosis factor-related protein 7
Entrez Gene ID	114905
mRNA Refseq	NM_001135170
Protein Refseq	NP_001128642
UniProt ID	Q9BXJ2