



CUEDC1 blocking peptide (CDBP5346)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CUEDC1 CUE domain containing 1 [Homo sapiens (human)]
Official Symbol	CUEDC1
Synonyms	CUEDC1; CUE domain containing 1; CUE domain-containing protein 1
Entrez Gene ID	404093
mRNA Refseq	NM_001271875
Protein Refseq	NP_001258804
UniProt ID	Q9NWM3