



CUEDC2 blocking peptide (CDBP5347)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CUEDC2 CUE domain containing 2 [Homo sapiens (human)]
Official Symbol	CUEDC2
Synonyms	CUEDC2; CUE domain containing 2; C10orf66; bA18I14.5; CUE domain-containing protein 2
Entrez Gene ID	79004
mRNA Refseq	NM_024040
Protein Refseq	NP_076945
UniProt ID	Q9H467
Function	protein binding