



CXXC4 blocking peptide (CDBP5353)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CXXC4 CXXC finger protein 4 [Homo sapiens (human)]
Official Symbol	CXXC4
Synonyms	CXXC4; CXXC finger protein 4; IDAX; CXXC-type zinc finger protein 4; CXXC finger 4; inhibition of the Dvl and axin complex protein; Dvl-binding protein IDAX (inhibition of the Dvl and Axin complex)
Entrez Gene ID	80319
mRNA Refseq	NM_025212
Protein Refseq	NP_079488
Pathway	Disease; RNF mutants show enhanced WNT signaling and proliferation; Signal Transduction; Signaling by WNT in cancer; Signaling by Wnt; TCF dependent signaling in response to WNT;

Wnt Signaling Pathway NetPath; Wnt signaling pathway

Function

DNA binding; PDZ domain binding; zinc ion binding
