



HSPC206 blocking peptide (CDBP5354)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	HSPC206 hypothetical protein HSPC206 [Homo sapiens (human)]
Official Symbol	HSPC206
Synonyms	HSPC206; hypothetical protein HSPC206
Entrez Gene ID	51525