



DACT2 blocking peptide (CDBP5355)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DACT2 dishevelled-binding antagonist of beta-catenin 2 [Homo sapiens (human)]
Official Symbol	DACT2
Synonyms	DACT2; dishevelled-binding antagonist of beta-catenin 2; DPR2; DAPPER2; C6orf116; bA503C24.7; dapper homolog 2; dapper antagonist of catenin 2; dapper homolog 2, antagonist of beta-catenin; dapper, antagonist of beta-catenin, homolog 2
Entrez Gene ID	168002
mRNA Refseq	NM_001286350
Protein Refseq	NP_001273279
UniProt ID	Q5SW24

Pathway	TGF-beta receptor signaling
Function	beta-catenin binding; delta-catenin binding; protein kinase A binding; protein kinase C binding; transcription factor binding