



DACT3 blocking peptide (CDBP5356)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DACT3 dishevelled-binding antagonist of beta-catenin 3 [Homo sapiens (human)]
Official Symbol	DACT3
Synonyms	DACT3; dishevelled-binding antagonist of beta-catenin 3; RRR1; DAPPER3; dapper homolog 3; arginine-rich region 1 protein; antagonist of beta-catenin Dapper homolog 3; dapper, antagonist of beta-catenin, homolog 3
Entrez Gene ID	147906
mRNA Refseq	NM_001301046
Protein Refseq	NP_001287975
UniProt ID	Q96B18

Function

delta-catenin binding; identical protein binding; protein kinase A binding; protein kinase C binding
