



## DBX1 blocking peptide (CDBP5361)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                               |
|-----------------|-----------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">DBX1 developing brain homeobox 1 [ Homo sapiens (human) ]</a>                     |
| Official Symbol | DBX1                                                                                          |
| Synonyms        | DBX1; developing brain homeobox 1; homeobox protein DBX1; developing brain homeobox protein 1 |
| Entrez Gene ID  | <a href="#">120237</a>                                                                        |
| mRNA Refseq     | <a href="#">NM_001029865</a>                                                                  |
| Protein Refseq  | <a href="#">NP_001025036</a>                                                                  |
| UniProt ID      | A6NMT0                                                                                        |
| Function        | sequence-specific DNA binding transcription factor activity; transcription regulatory region  |

