



Dbx2 blocking peptide (CDBP5362)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Dbx2 developing brain homeobox 2 [Mus musculus (house mouse)]
Official Symbol	Dbx2
Synonyms	Dbx2; developing brain homeobox 2; Gm1229; AI846217; 9630005K15; 9430056A22Rik
Entrez Gene ID	223843
mRNA Refseq	NM_207533
Protein Refseq	NP_997416
Function	DNA binding