



DCLK3 blocking peptide (CDBP5364)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DCLK3 doublecortin-like kinase 3 [Homo sapiens (human)]
Official Symbol	DCLK3
Synonyms	DCLK3; doublecortin-like kinase 3; CLR; DCK3; DCDC3C; DCAMKL3; serine/threonine-protein kinase DCLK3; CLICK-I,II-related; doublecortin and CaM kinase-like 3; doublecortin-like and CAM kinase-like 3; doublecortin domain-containing protein 3C
Entrez Gene ID	85443
mRNA Refseq	NM_033403
Protein Refseq	NP_208382
UniProt ID	Q9C098

Function

ATP binding; protein serine/threonine kinase activity
